

TEACHING THE
GENOME
GENERATION

PREP FOR SEQUENCING PROTOCOL

PREREQUISITES & GOALS

PREREQUISITES

Prior to implementing this lab, you should understand:

- The central dogma of how DNA bases code for mRNA and then for proteins
- How DNA samples were collected and prepared for PCR
- The steps that occur during the process of polymerase change reaction (PCR), as well as the reagents used
- The variants of ACTN3, CYP2C19, and/or TAS2R38
- The connection between the bands on the gel electrophoresis and the sequential gene variants
- The purpose of using ExoSAP-IT
- The purpose of the PREP FOR SEQUENCING PROTOCOL is to prepare PCR-amplified DNA samples for Sanger sequencing

LEARNING GOALS

1. Prepare amplified DNA products for Sanger sequencing.

NOTES

The **ExoSAP** system uses two enzymes to prepare amplified PCR products for sequencing:

Exonuclease I (Exo) digests residual single-stranded DNA primers and any extraneous single-stranded DNA fragments produced in from the PCR process

Shrimp Alkaline Phosphatase (SAP) dephosphorylates remaining dNTPs (free nucleotides) from the PCR product so they do not interfere with the sequencing reaction

MATERIALS

REQUIRED LAB MATERIALS

Ice bath or crushed ice

Markers for labeling

Amplified DNA samples from the PCR PROTOCOL

Micropipettors & tips
(size P20)

0.2 mL tubes in strips of 4

1.5 mL tubes

Tube holders/racks

Mini-microcentrifuge

Thermal cycler

Molecular biology grade water

(F)orward sequencing primer
(4 μ M)

(R)everse sequencing primer
(4 μ M)

Invitrogen ExoSAP-IT enzyme
(on ice)

WORKSTATION NEEDS

These materials should be at each workstation

Micropipettors and tips

0.2 mL tubes in strips

1.5 mL tubes

Molecular biology grade water

(F)orward and (R)everse primers

Tube holders and markers
for labeling

Crushed ice/ice bath

Exo-SAP IT enzyme (on ice)

Amplified DNA samples

PROCEDURE

□ STEP 1

Obtain 0.2 mL strip tubes and label them with the DNA sample numbers/letters.

NOTE: Do not use the negative control sample in this procedure, only samples that demonstrated positive amplification during gel electrophoresis should be processed.


□ STEP 2

Using the P20 micropipettor, transfer 5 μ L of the amplified DNA samples from PROTOCOL 2 to fresh tubes.

□ STEP 3

Using the P20 micropipettor, add 2 μ L of ExoSAP-IT to the new tubes.

NOTE: Enzymes must be kept on ice at all times.



NOTES

A large grid of dots for taking notes.

□ STEP 4

Tightly cap the tubes and flick the tube gently to mix.

□ STEP 5

Spin the tubes briefly in the mini-microcentrifuge to collect the solution.

BREAK POINT IF NEEDED.

Expected result is to have one tube per DNA sample with 7 μL of reaction solution.

□ STEP 6

The thermal cycler should be programmed with the digestion protocol.

EXOSAP: Digests free nucleotides and primers

Cycling conditions

- | | | |
|------------------------|-------|---------|
| 1. Digestion | 37° C | 15 min. |
| 2. Protein degradation | 85° C | 15 sec. |
| 3. Final hold | 4° C | forever |

□ STEP 7

Consult your teacher on proper use of the thermal cycler provided.

BREAK POINT - the reaction will proceed for 16 minutes.

NOTES

STEP 8

Remove the samples after the protocol is complete, stop the program and turn the machine off.

BREAK POINT IF NEEDED.

Expected result is to have one tube per DNA sample with 7 μL of reaction solution. Nothing should look different about the solution after the reaction.

STEP 9 - CRITICAL STEP

Obtain two new 0.2 mL strip tubes per DNA sample to be sequenced.

STEP 10 - CRITICAL STEP

Label one tube as the forward sequencing reaction and one as the reverse sequencing reaction for each sample.


EXAMPLE: 4-F would be the label for DNA sample #4 to be sequenced with the Forward primer.

STEP 11

Using the P20 micropipette, add the following to the FORWARD reaction tube:

- 10 μL of molecular biology grade water
- 3 μL of ExoSAP-IT/DNA mix that has been incubated
- 2 μL of (F)orward sequencing primer

NOTES



□ STEP 12

Using the P20 micropipettor, add the following to the REVERSE reaction tube:

- 10 μ L of molecular biology grade water
- 3 μ L of ExoSAP-IT/DNA mix that has been incubated
- 2 μ L of (R)everse sequencing primer

□ STEP 13

Tightly cap the tubes and flick the tube gently to mix.

□ STEP 14

Spin the tubes briefly in the mini-microcentrifuge to collect the solution.

Expected result is to have two tubes per DNA sample with 15 μ L of reaction solution.

□ STEP 15

Follow the final sample labeling, shipping and packaging directions dictated by the sequencing company.

The DNA sample is now ready for sequencing.

NOTES

