

Updated 02/09/2024

# TEACHING THE GENOME GENERATION™

## *DNA EXTRACTION PROTOCOL (LONG)*

### BEFORE YOU BEGIN

Have you discussed *informed consent* with your students? This key feature of the TtGG curriculum is vitally important in their decision whether to provide a saliva sample for downstream processing. We strongly encourage you go over this important subject prior to beginning the laboratory experiments.

# PREREQUISITES & GOALS

## STUDENT PREREQUISITES

Prior to implementing this lab, students should understand:

- The central dogma of how DNA bases code for mRNA and then for proteins
- The purpose of the DNA EXTRACTION PROTOCOL is to extract human DNA and make the sample ready to amplify in the PCR PROTOCOL
- How the saliva kit is extracting DNA
- Units of measurement ( $\mu\text{L}$ )

## STUDENT LEARNING GOALS

1. Complete lab procedures necessary to collect DNA samples.
2. Perform proper micropipetting technique.
3. Identify ethical issues with DNA sample collection.
4. Extract and see a pellet of DNA.
5. Use lab equipment, including the centrifuge and the vortex.

## ETHICAL ISSUES

This protocol uses saliva and cheek cells as a source for extracting purified human DNA. All experiments in the course are demonstrations; none of the genotyping performed on the human samples are in any way diagnostic. For a number of ethical reasons, it is very important to allow the DNA collection stage to be 100% voluntary. There are personal, cultural, religious, and privacy based reasons why students may not want to participate.

Although most students will want to know their own personal genotype or DNA sequence, don't give in to their pleadings. It is imperative that the samples collected are not labeled by name, number or category of any kind. The goal is to keep samples anonymous and not be able to match sample to person. At the end of this protocol, collect the unlabeled DNA tubes and put generic labels (1, 2, 3, etc. or A, B, C, etc.) on the tubes prior to starting subsequent procedures.

# CURRICULUM INTEGRATION

Use the planning notes space provided to reflect on how this protocol will be integrated into your classroom. You'll find every course is different, and you may need to make changes in your preparation or setup depending on which course you are teaching.

Course name:

1. What prior knowledge do the students need?

2. How much time will this lesson take?

3. What materials do I need to prepare in advance?

4. Will the students work independently, in pairs, or in small groups?

5. What might be challenge points for students during this lesson?

# MATERIALS

## REQUIRED LAB MATERIALS

Ice bath or crushed ice

Markers for labeling

Refrigerator

## PROVIDED BY JAX

Provided for TtGG-trained teachers, contact [ttgg@jax.org](mailto:ttgg@jax.org).

DNA Genotek Oragene  
DNA collection kits

Oragene-DNA "OG-L2P"  
purifier solution

Micropipettes & tips  
(sizes P1000, P200 & P20)

1.5 mL tubes

Tube holders/racks

Centrifuge

Ethanol (100% and 70%)

Molecular biology grade water

Heat Block

Vortex

## WORKSTATION NEEDS

*Distribute these materials to each workstation.*

DNA collection kit and purifier solution

Micropipettes and tips

1.5 mL tubes

Ethanol (100% and 70%)

Molecular biology grade water

Tube holders

Markers for labeling

Crushed ice/ice bath

# PROTOCOL STRUCTURE

**STEPS 1** 8 minute video on proper pipetting technique *(critical for good results)*

Break point

**STEPS 2-8** 20 minutes

Break point: samples can be stored at 4 °C for several weeks

**STEP 9** 90 minutes minimum *(incubation period – students need not be present)*

Break point: samples can be stored at 4 °C for several weeks

**STEPS 10-30** 60-70 minutes

Break point: sample can be dried overnight

**STEPS 31-32** 5 minutes

# PROCEDURE

## STEP 1

Watch: Using a Micropipette — University of Leicester  
[www.youtube.com/watch?v=uEy\\_NGDfo\\_8&sns=em](http://www.youtube.com/watch?v=uEy_NGDfo_8&sns=em)  
and/or perform the MICROPIPETTING EXERCISE.

NOTE: Micropipetting skill is a critical component of the lab. If students do not do this part well, the results of all subsequent activities may be compromised.

## BREAK POINT IF NEEDED

## STEP 2

Obtain an Oragene-DNA saliva kit.

## STEP 3

Each student needs to produce enough saliva to fill the sample tube up to the black line with fluid — not bubbles. (Yes, students spit into the tube!)

## STEP 4

Once the saliva is in the tube, close the cap. This releases the required buffer.

*WHY: This proprietary buffer lyses cells and stabilizes the DNA in solution.*

## STEP 5


After the buffer is released, the funnel cap should be removed and discarded and the small circular cap should be put on the tube.

## STEP 6

Mix the saliva in the sample tube by inversion and gentle shaking for a few seconds.

NOTE: DNA in Oragene-DNA buffer is stable at room temperature and can be kept as backup until all procedures are complete.

## PLANNING NOTES



## □ STEP 7

Obtain a 1.5 mL microcentrifuge tube. Label the tube with the student's initials.

## □ STEP 8

Use the P1000 micropipette to transfer 500  $\mu$ L of the mixed Oragene-DNA/saliva sample into the tube.

NOTES:

- If there is mucus-like material in the collection tube, try to avoid sucking up the viscous mucus component.
- There may be color in the solution. For instance, if you collect saliva while drinking coffee the solution will be brown. It is best to collect saliva when not eating or drinking.

## BREAK POINT IF NEEDED

*Samples can be stored at room temperature for a few hours, 4 °C (refrigerator) for several weeks, or frozen for long term storage.*

## □ STEP 9

Incubate the sample tubes at 50 °C in the heat block for a minimum of 90 minutes.

*WHY: This heating step will denature the proteins in the solution.*

NOTE: Samples can incubate overnight at 50 °C if necessary, but this is not recommended.

## BREAK POINT IF NEEDED

*Samples can be stored at 4 °C (refrigerator) for several weeks.*

## PLANNING NOTES

A large grid of dotted lines for planning notes, consisting of 20 columns and 30 rows of small grey dots.

## □ STEP 10

Return each student's respective sample.

## □ STEP 11

Using a P20 micropipette, add 20  $\mu$ L Oragene-DNA "OG-L2P" purifier to each sample and mix by vortexing for a few seconds.

NOTE: OG-L2P purifier tubes should be labeled with students' initials.

## □ STEP 12

Incubate on ice for 10 minutes.

*While waiting, read ahead and prepare tubes for Step 16.*

*WHY: The purifier reagent and cold temperature will precipitate cell debris (proteins, lipids, etc.) from the solution. DNA will remain in solution.*

## □ STEP 13

Place tubes in centrifuge with the hinge of the tube along the outer edge of the rotor (see photo). The tubes need to be balanced in the centrifuge. Tubes with water can be used to balance, if needed.

NOTE: Make sure tubes are labeled before they enter the centrifuge so the students can retrieve their own samples after the spin.



## PLANNING NOTES

A large grid of small dots for planning notes.

## □ STEP 14

Centrifuge at room temperature for 5 minutes at 13,000 rpm (15,000 x g).

*WHY: Centrifugation pellets (or physically separates) the precipitated cell debris from the DNA in solution.*

## □ STEP 15

Remove tubes from centrifuge and look for a pellet on the bottom or side of the tube, on the side of the hinge.

NOTES:

- a. If a pellet doesn't form or is just a viscous mess, re-centrifuge for up to 15 minutes.
- a. The pellet contains insoluble material and impurities. The DNA is in the soluble supernatant (liquid phase) and must be transferred into a fresh tube.


## □ STEP 16

Obtain and label a fresh 1.5 mL tube.

## □ STEP 17

Using the P1000 micropipette, carefully transfer 400  $\mu$ L of the clear liquid supernatant into the tube.

**THE CLEAR SUPERNATANT CONTAINS THE DNA.**



## PLANNING NOTES

A large grid of dots for planning notes.



## NOTES:

- For a more precise supernatant transfer, the P200 can be used with two draws.
- Do not disturb the pellet. However, if the pellet is disturbed, re-centrifuge and transfer again.
- Discard the tube with the pellet in the trash.

## BREAK POINT IF NEEDED

*Samples can be stored at 4 °C for 1-2 days.*

 STEP 18

Add 400  $\mu\text{L}$  of room temperature 100% ethanol to the 400  $\mu\text{L}$  of the supernatant.

*WHY: DNA is hydrophilic, so adding ethanol will allow DNA to precipitate out of the solution.*

NOTE: The volume of ethanol added must be equal to the volume of supernatant, so if you have less than 400  $\mu\text{L}$  of supernatant, add less ethanol.

 STEP 19

Mix the tube gently by inversion ~10 times.

NOTE: The DNA may precipitate out and appear as a clot of fibers or cloudiness.

 STEP 20

Let stand for ~5 minutes at room temperature to allow the DNA to fully precipitate.

 STEP 21

Centrifuge at room temperature for 2 minutes at 13,000 rpm (15,000  $\times g$ ), ensuring the centrifuge is properly balanced (see STEP 13).

*WHY: The precipitated DNA is physically separated from the solution by centrifugation.*

## PLANNING NOTES



## □ STEP 22

Remove the tube and locate the position of the DNA pellet. It should be at the inside bottom portion of the tube below the hinge. (In circle area of tube, below hinge, see photo.)



### NOTES:

- Now the pellet contains the DNA. The pellet is very small and will not always be visible, this is okay. Proceed to next steps.
- If the pellet is disturbed, re-centrifuge.
- Now the supernatant contains impurities (not the DNA) and will be discarded in a later step.

## □ STEP 23

Rotate the tube so the DNA pellet and hinge are on the upper wall of the tube (see diagram above). This will allow easier access to extract the supernatant.

## □ STEP 24

Dial the P1000 micropipette to 800  $\mu$ L and carefully remove and discard the supernatant. **THE PELLETT CONTAINS YOUR DNA.**

### NOTES:

- It is not necessary to remove all of the supernatant in this step.
- Do not disturb the DNA pellet with the pipette tip.
- The P200 micropipette can be used with multiple draws to reduce the risk of pellet disturbance

## PLANNING NOTES

A large grid of small dots for planning notes, consisting of 20 columns and 30 rows of light blue dots on a white background.

and allow for a more precise extraction of the supernatant. Dial the pipettor to 200  $\mu\text{L}$ .

## STEP 25

Using the P1000 micropipette, add 250  $\mu\text{L}$  of 70% ethanol to wash the DNA pellet.

*WHY: The DNA pellet needs to be washed with ethanol again to ensure all impurities are removed.*

## STEP 26

Centrifuge at room temperature for 2 minutes at 13,000 rpm (15,000  $\times$  g), ensuring the centrifuge is properly balanced (see STEP 13).

## STEP 27

Remove the tube and locate the DNA pellet. The pellet is very small and will not always be visible, this is okay. Proceed to next steps.

## STEP 28

Using the P200 micropipette, remove all supernatant without disturbing the DNA pellet.

NOTE: To remove all of the supernatant, it will take two draws with the micropipette.

## STEP 29

Carefully turn tube upside down over a paper towel to allow excess drops of ethanol to run out.

## STEP 30

Keep the tube open, place it on its side and store at room temperature to allow ethanol to evaporate completely for 10-30 minutes. Tubes can be left open overnight as well.

*WHY: All ethanol needs to evaporate or the DNA will not go back into solution fully; residual ethanol may hinder downstream reactions.*

## PLANNING NOTES



## BREAK POINT IF NEEDED

*Samples can be dried overnight if needed and stored at 4 °C (refrigerator) for several weeks.*

### STEP 31

After the pellet is dry, use the P200 micropipette to add 100  $\mu$ L of molecular biology grade water to dissolve the DNA pellet.

*WHY: Molecular Biology-grade water is very pure and will not inhibit downstream reactions, unlike other mineralized or ionized water.*

### STEP 32

Wait 5 minutes and vortex at medium speed for 10 seconds, or until the pellet has dissolved. The expected concentration of the fully hydrated DNA is 20-200 ng/ $\mu$ L.


Expected result is to have one tube with 100  $\mu$ L of pure, hydrated DNA per saliva sample.

The DNA sample can be stored at -20 °C (freezer) for up to 5 years.

Collect student samples and remove any labels on the tubes by wiping off the marker with left-over ethanol.

Re-label tubes with arbitrary letters or numbers to preserve anonymity.

## PLANNING NOTES



## The DNA sample is now ready for POLYMERASE CHAIN REACTION (PCR) PROTOCOL

### Sources of Potential Error:

The most common error for DNA EXTRACTION PROTOCOL (LONG) is incorrect micropipetting or accidental loss of the pellet.

### Clean up:


Unused DNA sample collection kits should be saved and returned to JAX.

Any remaining sample from STEP 6 can be stored at room temperature, or in the freezer for longer term storage, to repeat this assay another time, if desired.

Discard all used tubes except the hydrated DNA sample.

All used tips and tubes can be disposed of in the trash.

## PLANNING NOTES

A large grid of small dots, intended for planning notes. The grid consists of 20 columns and 30 rows of dots, providing a structured space for students to write down their experimental plans.

## NEED HELP?

Email the experts: [ttgg@jax.org](mailto:ttgg@jax.org)